

AUTOMATIC EXTERNAL DEFIBRILLATOR



A
E
D

Accurate and Reliable ECG Analysis
Easy to Use
Motion Detector
Pace Maker Detection

Self Test Function
AHA/ERC Guidelines for CPR and Defibrillation
5 Year Warranty



Turkish First Aid Association supports "Life-Point"

CE 1984



**AUTOMATIC
EXTERNAL
DEFIBRILLATOR**

Life-Point Pro AED Technical Specifications

Operating Mode	Automatic
Charging Time at max Joule	typically < 10 saniye
ECG Analysis Period	typically < 10 saniye
Wave Form	Truncated Biphasic Wave Form
Energy Steps	Adult :150J - 150J - 200J / Child : 50J - 50J - 50J
Pace Maker Detection	Yes
Motion Detector	Yes
Impedance Measurement	25 – 100 ohm
Maximum Discharging	typically 200 shocks (at 20 °C with full battery)
Derivation	D II
ECG Bandwidth	2 – 25 Hz
Battery Type	LiMnO ₂ (15V 3000mA)
Self Test Schedule	Weekly - Monthly - Each Start
Self Test Scope	Battery - Internal Electronics - Software
Dimensions	9,5cm x 27cm x 24cm
Weight (battery included)	1730 gr
Class (93/42/EEC)	II B
Operating Conditions	0°C - 50°C ve % 10 – % 95 RH
Storage Conditions	Between -20°C and 60°C
Data Documentation	ECG Signal/ Event / Voice Record
Data Transfer	Metsis AED Software (USB/SD Kart)
Minimum Computer Requirement	Windows XP – Windows Vista – Windows 7 Pentium or higher VGA Monitor or higher CD-ROM Driver USB Port 2 Gb Free Disk Space

www.metsismedikal.com



METSis Medikal Elektronik Ltd. Şti.

Adres : İvedik Organize Sanayi Bölgesi
Soğutmaçılar Sitesi 1444.Sokak No:11
06370 Ostim - Ankara - TURKEY

Tel : +90 312 395 77 70

Fax : +90 312 395 77 72

e-mail : info@metsismedikal.com